



Little Step-by-Step

Product Information	
Item Number	100-02500 : Includes Red, Yellow, Blue and Green
Cost	\$146
Memory	120 seconds variable sequential message location
Applications	Communication
Size / weight	4"x4"x3" (10.16cmx10.16cmx7.62cm); 8.5 oz (241g)
Additional Information	Now includes changeable color switch top (set of 4)
Approvals	CE, RoHS
Medicare version/code	E2500
Device Groups	
Device type	Digital recorded voice output
Input/Selection Method	
Selection method	Direct select with external Input jack
Feedback	Message (Entire message play regardless of repetitive hits)
Activation Surface	Tactile Switch with Activation Pressure 2-3 oz (57-85 grams)
Types of scanning	None
Size of message Locations	2.5" (6.35 cm)
Number of message locations	1 (Variable sequential messages)
Numbers of levels	None
Speech output	
Speech output	Digital recorded voice output
Advance Features	
	Repeat feature; playing a message with an external switch; entire message play regardless of repetitive hits; message can activate toy/appliance
Accessories	
Accessories	Includes snap switch cap and toy/appliance cable and easy mounting base
	Optional: AbleNet Belt Mount strap, Small Triangula mounting Plate
Power storage	
Battery information	9 volt battery (not included); 4.25 hours of continues talk time; Sleep mode
Hardware info	
Speaker / microphone / headphone	Built-in Speaker and microphone
Controls	
Programming	Record
Volume control	Dial
Volume Lock Out	No
Power on/off	Dial (Part of Volumn Control)
Warranty information	
Warranty information	1 year