

[-- Back to Switches, Computer Access, Controls](#)

## Switches

Small, soft, sensitive, strong, ribbon, lever, scatic, dual, multiple, & adjustable pressure switch with buddy button



Small Switches



Soft Switches




AbleNet Access Switches



Sensitive Switches




Buddy Buttons



Ribbon Switches




Lever Switches




Strong Switches



Dual Switches



Multiple Switches



Adjustable Pressure Switch



Accessories/Cables/Kits

Using PayPal a problem? Fax your purchase order to 830-608-0882.  
 For discount on shipping multiple items call 877-608-9812.  
[Link to calculate shipping for P.O.'s](#)  
[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.  
 Prices subject to change without notice.

Turning Point Therapy & Technology,  
 Inc.  
 PO Box 310945  
 New Braunfels, TX 78131-0945  
 Member of Assistive Technology Industry  
 Association

**Contact Us**  
 Toll Free: 877-608-9812  
 Phone/Fax: 830-608-9812/830-608-0882  
 Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
 Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
 Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)

Since 11/24/98



[← Back to Switches](#)

## Small Switches

We have a wide range of single switches used for people with a wide range of capabilities, all with an 1/8 in mono plug and a 5' cable, including small trigger, cap, cup, heavy duty cap, and minicap switch

Item	Description	Photo	Price
Switch <a href="#">TA-58060</a>	<b>Cap Switch and Mount/Connection</b> <ul style="list-style-type: none"> <li>• Activated by pressing down on the top surface.</li> <li>• Has a non-slip base for mounting flat.</li> <li>• Tactile and auditory feedback.</li> </ul> <a href="#">Link to Mounting Options</a>	Switch 	Switch <b>\$63.00</b>
Connection <a href="#">TA-66110</a>		Connection 	Connection <b>\$27.00</b>
Switch <a href="#">TA-58070</a>	<b>Cup Switch, Buddy Button and Mount/Connection</b> <ul style="list-style-type: none"> <li>• Same surface area as the cap, but more durable since the switch is encased in a plastic cup.</li> <li>• Comes with colored stickers and Velcro for easy mounting.</li> <li>• Tactile and auditory feedback.</li> </ul> <a href="#">Link to Mounting Options</a>	Switch 	Switch <b>\$74.00</b>
Connection <a href="#">TA-66120</a>		Connection 	Connection <b>\$37.00</b>
<a href="#">TA-58830</a>	<b>Heavy Duty Cap Switch</b> <ul style="list-style-type: none"> <li>• A heavy duty plastic makes this switch extremely rugged and superior in quality.</li> <li>• Activated by pressing down anywhere on the top surface.</li> <li>• Tactile and auditory feedback.</li> </ul> <a href="#">Link to Mounting Options</a>		<b>\$158.00</b>

Switch  
TA-58080  
Connection  
TA-66110

### Mini Cup Switch and Mount/Connection

- Very small switch for people with good targeting ability; activated by pressing on surface with finger, chin, etc.
- Comes with colored stickers.
- Tactile and auditory feedback.

[Link to Mounting Options](#)



Switch  
**\$63.00**  
Connection  
**\$27.00**

Connection



**gray**  
TA-55000

**red**  
TA-55100

**green**  
TA-55200

**orange**  
TA-55400

**blue**  
TA-55600

### Trigger Switch

- Designed to meet the rugged needs of wheelchair users.
- Waterproof and dust proof, ideal for rugged users in any environment.
- The design accommodates mounting on a flat surface with all mounting solutions included in the price.
- Comes in 5 colors

[Link to Mounting Options](#)



Switch  
**\$64.00**

Using PayPal a problem? Fax your purchase order to 830-608-0882.  
 For discount on shipping multiple items call 877-608-9812.  
[Link to calculate shipping for P.O.'s](#)  
 Shipping outside the continental United States requires [email contact](#) prior to purchase.  
 Prices subject to change without notice.

Turning Point Therapy & Technology,  
 Inc.  
 PO Box 310945  
 New Braunfels, TX 78131-0945

Member of Assistive Technology Industry  
 Association

#### Contact Us

Toll Free: 877-608-9812  
 Phone//Fax: 830-608-9812//830-608-0882  
 Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
 Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
 Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)

Since 11/24/98







[← Back to Switches](#)

## Soft Switches

We have a wide range of single switches used for people with a wide range of capabilities, all with an 1/8 in mono plug and a 5' cable

Item	Description	Photo	Price
TA-58750	<p>Pillow Switch</p> <p>Activated by pressing the top foam surface. Smooth, soft surface makes this a suitable switch for head or cheek activation. Safety pin and velcro come on the bottom for easy fastening to a pillow or wheelchair cushion. Tactile and auditory feedback.</p>		\$105.00
<p>Red TA-58200</p> <p>Green TA-58210</p>	<p>Soft Switches</p> <p>An excellent switch for people with gross motor skills. The switch is enclosed in foam and covered with a removable, washable velvet bag that comes in 2 colors: Tactile and auditory feedback.</p>		\$84.00

Using PayPal a problem? Fax your purchase order to 830-608-0882.

For discount on shipping multiple items call 877-608-9812.

[Link to calculate shipping for P.O.'s](#)

[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.

Prices subject to change without notice.

Turning Point Therapy & Technology,  
Inc.  
PO Box 310945  
New Braunfels, TX 78131-0945

Member of Assistive Technology Industry  
Association

### Contact Us

Toll Free: 877-608-9812  
 Phone//Fax: 830-608-9812//830-608-0882  
 Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
 Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
 Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)

Since 11/24/98





[← Back to Switches](#)

## AbleNet Access Switches



Item	Description	Photo Click Photo to Enlarge	Price
AN-1-ACAL AirLink Cordless Switch	<p>2 switches in one! Cordless infrared access to the PowerLink 3 control unit allows students to activate appliances up to 20 feet away. No cords means no hassle during turn-taking activities and easy access for all students even in small classrooms. Attach AirLink's cord appliance cable (included) to create corded switch access. AirLink features a 2 1/2-inch activation surface on an angled base, and is recommended for students who can access a smaller target area.</p> <p>Features:</p> <ul style="list-style-type: none"> <li>• External switch jack for specialty switch operation</li> <li>• Requires 9-volt battery (not included)</li> <li>• Switch access cable included</li> </ul>		<b>\$89.00</b>
Red AN-1-BR Black AN-1-BBK Blue AN-1-BBL Green AN-1-BG Yellow AN-1-BY	<p><b>Big Red Switch</b></p> <p>The Big Red switch features a large, 5-inch activation surface, and is perfect for students who require a larger surface due to visual and motor abilities.</p>		<b>\$45.00</b>
Red AN-100-30101 Blue AN-100-30103 Green AN-100-30105 Yellow AN-100-30106	<p><b>Jelly Beamer Switch</b> Jelly Bean Switch... Without the Wires</p> <ul style="list-style-type: none"> <li>• Rugged, reliable and easy to use, designed and built with the same hardwearing features you've come to expect from all AbleNet switches.</li> <li>• Uses the same technology found in wireless headsets for mobile phones and MP3 players, &amp; PDAs.</li> <li>• Simple Set-up. Immediate Action: just 2 components: the switch that houses a transmitter &amp; the Jelly Beamer receiver that can be plugged into any switch-adapted device. Each switch is uniquely paired with its receiver.</li> <li>• Advantages: Easily be passed between students, Fosters turn-taking and sharing, Reduce cord clutter in the work area</li> <li>• Quantity pricing makes buying multiple Jelly Beamer's affordable:</li> </ul> <p>1 - 2 switches = \$129.00 each 3 or more switches = \$99.00 each</p>		<p><b>\$135.00</b></p> <p>Call 877-608-9812 for Quantity Discount</p>
<a href="#">Details</a>			

Red  
AN-1-JBR

### Jelly Bean Switch

\$45.00

Considered a standard in the industry, AbleNet's Jelly Bean switch features a 2 1/2-inch activation surface and is recommended for students who can access a smaller target area.



Black  
AN-1-JBBK

Blue  
AN-1-JBBL

Green  
AN-1-JBG

Yellow  
AN-1-JBY

Red  
AN-1-SPR

### Specs Switch

\$49.00

The Specs switch is the most versatile and durable small switch available. With a 1 3/8-inch activation surface, this switch is small enough to be worn, but is most often used as a mounted switch. The Specs switch includes 3 bases: a standard flange base, a space saving flush base and a strap base for mounting around items. 24-inch strap included.



Black  
AN-1-SPBK

Blue  
AN-1-SPBL

Green  
AN-1-SPG

Yellow  
AN-1-SPY

Using PayPal a problem? Fax your purchase order to 830-608-0882.  
For discount on shipping multiple items call 877-608-9812.  
[Link to calculate shipping for P.O.'s](#)  
[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.  
Prices subject to change without notice.

Turning Point Therapy & Technology,  
Inc.  
PO Box 310945  
New Braunfels, TX 78131-0945

Member of Assistive Technology Industry  
Association

#### Contact Us

Toll Free: 877-608-9812  
Phone//Fax: 830-608-9812//830-608-0882  
Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
Website: [webmaster@turningpointtechnology.com](http://webmaster@turningpointtechnology.com)





Since 11/24/98



[← Back to Switches](#)

## Sensitive Switches

We have a wide range of single switches used for people with a wide range of capabilities, all with an 1/8 in mono plug and a 5' cable, including sensitive leaf, plate, tip, and microlight switch

Item	Description	Photo	Price
TA-58600	<p>Leaf Switch</p> <ul style="list-style-type: none"> <li>• Activated by pressing lightly against the leaf in one direction.</li> <li>• Comes with a removable 3" round foam pad.</li> <li>• Tactile and auditory feedback.</li> <li>• Requires mounting.</li> </ul>		<b>\$147.00</b>
<p>Switch TA-58500</p> <p>Connection TA-66100</p>	<p>Microlight Switch and Mount/Connection</p> <ul style="list-style-type: none"> <li>• Activated by pressing down on the top surface with a feather light touch!</li> <li>• Threaded insert on the bottom makes mounting easy with a standard "4-40" screw (available at electronics stores). Tactile and auditory feedback.</li> </ul> <p><a href="#">Link to Mounting Options</a></p>	<p>Switch</p>  <p>Connection</p> 	<p>Switch <b>\$69.00</b></p> <p>Connection <b>\$37.00</b></p>
<p>Switch <b>red</b> TA-58300</p> <p><b>orange</b> TA-58340</p> <p><b>yellow</b> TA-58320</p> <p><b>green</b> TA-58310</p> <p><b>blue</b> TA-58330</p> <p>Connection TA-66130</p>	<p>Plate Switch and Mount/Connection</p> <ul style="list-style-type: none"> <li>• A very thin membrane switch activated with a light touch on the colored circle.</li> <li>• Surface is sealed to prevent moisture damage.</li> <li>• No tactile/auditory feedback.</li> <li>• Comes in 5 bright colors</li> </ul> <p><a href="#">Link to Mounting Options</a></p>	<p>Switch</p>  <p>Connection</p> 	<p>Switch <b>\$69.00</b></p> <p>Connection <b>\$37.00</b></p>

TA-58400

## Tip Switch

- A mercury switch that is activated by tilting 5 degrees from horizontal position.
- Great for mounting on a hand or head band - requires a very small range of motion.
- No tactile/auditory feedback.



\$69.00

Using PayPal a problem? Fax your purchase order to 830-608-0882.

For discount on shipping multiple items call 877-608-9812.

[Link to calculate shipping for P.O.'s](#)

[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.

Prices subject to change without notice.

Turning Point Therapy & Technology,  
Inc.  
PO Box 310945  
New Braunfels, TX 78131-0945

Member of Assistive Technology Industry  
Association

### Contact Us













Toll Free: 877-608-9812  
Phone//Fax: 830-608-9812//830-608-0882  
Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)

Since 11/24/98





[← Back to Switches](#)

## Buddy Buttons

Item	Description	Photo Click to enlarge	Price
 Black TA-57000	<b>Buddy Buttons</b> <ul style="list-style-type: none"> <li>The reliable, durable, colorful single switch, activated by pressing anywhere on the top surface.</li> <li>Auditory and tactile feedback.</li> <li>Available in standard colors</li> </ul>		Buddy Button <b>\$53.00</b>
 Blue TA-57600			
 Green TA-57200	<b>Specifications</b> Diameter 2.5/6.3 cm Height 0.7 in/1.8 cm Weight 2 oz/60 g		Undercover Replacement Stickers (4) <b>\$10.00</b>
 Hands On TA-57740			
 Happy TA-57720	<b>Activating Force</b> Outside edges 3.5 oz/100 g Center 7 oz/200 g Travel before activation 0.05 in/12 cm		
 Pink TA-57400			
 Purple TA-57300	<b>Current</b> 30VDC max 1A		
 Red TA-57100			
 Undercover TA-57800			
 Yellow TA-57500			
 White TA-57050			

## Mounts/Connections

TA-66120	Buddy Button Connection		<b>\$37.00</b>
TA-66130	Big Buddy Connection		<b>\$37.00</b>

Using PayPal a problem? Fax your purchase order to 830-608-0882.

For discount on shipping multiple items call 877-608-9812.

[Link to calculate shipping for P.O.'s](#)

[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.

Prices subject to change without notice.

Turning Point Therapy & Technology,  
Inc.  
PO Box 310945  
New Braunfels, TX 78131-0945

Member of Assistive Technology Industry  
Association

### Contact Us

Toll Free: [877-608-9812](tel:877-608-9812)  
Phone//Fax: [830-608-9812](tel:830-608-9812)//[830-608-0882](tel:830-608-0882)  
Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)

Since 11/24/98





[← Back to Switches](#)

## Ribbon Switch

We have a wide range of single switches, including ribbon switch, used for people with a wide range of capabilities, all with an 1/8" mono plug and a 5' cable.

Item	Description	Photo	Price
Ribbon Switch TA-58700	<p><a href="#">Ribbon Switch &amp; Mounting Arm</a> <a href="#">More on Mounting Options</a></p> <ul style="list-style-type: none"> <li>• Activated by bending in either direction</li> <li>• Soft to touch</li> <li>• Moisture resistant</li> <li>• Waterproof</li> <li>• Tactile feed back</li> <li>• Consistent positive action</li> </ul> <p><a href="#">Specifications</a> Size: 13 x 1.5 x .5 mm / 5 x 0.5 x 0.25 in Force: 110 g / 4 oz Weight: 32 g / 1.1 oz</p>		Ribbon Switch <b>\$105.00</b>

Using PayPal a problem? Fax your purchase order to 830-608-0882.

For discount on shipping multiple items call 877-608-9812.

[Link to calculate shipping for P.O.'s](#)

[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.

Prices subject to change without notice.

Turning Point Therapy & Technology,  
Inc.  
PO Box 310945  
New Braunfels, TX 78131-0945

Member of Assistive Technology Industry  
Association

### Contact Us

Toll Free: [877-608-9812](tel:877-608-9812)  
Phone//Fax: [830-608-9812](tel:830-608-9812)/[830-608-0882](tel:830-608-0882)  
Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)



Since 11/24/98



[← Back to Switches](#)

## Lever Switches

A wide range of single switches, including lever switch, used for people with a wide range of capabilities, all with an 1/8" mono plug and a 5' cable

Item	Description	Photo	Price
Black TA-53100	<b>MiniLever Switch</b> <ul style="list-style-type: none"> <li>• Simple Concept</li> <li>• Robust Construction</li> <li>• Auditory Feedback</li> <li>• Tactile Feedback</li> <li>• Available in the illustrated colors</li> <li>• Specifications given below</li> </ul>		<b>\$74.00</b>
Red TA-53110			
Green TA-53120			
Yellow TA-53150	<b>Specifications</b> Diameter: 107 x 37 x 18 mm / 4 x 1.5 x 0.75 in Force: 40 g / 1.4 oz Weight: 45 g / 1.5 oz		
Blue TA-53160			
Black TA-53000	<b>MicroLever Switch</b> <ul style="list-style-type: none"> <li>• Simple Concept</li> <li>• Robust Construction</li> <li>• Auditory Feedback</li> <li>• Tactile Feedback</li> <li>• Available in the illustrated colors</li> <li>• Specifications given below</li> </ul>		<b>\$72.00</b>
Red TA-53010			
Green TA-53020			
Yellow TA-53050	<b>Specifications</b> Diameter: 40 x 36 x 20 mm / 1.6 x 1.5 x 0.7 in Force: 40 g / 1.4 oz Weight: 31 g / 1 oz		
Blue TA-53060			

Using PayPal a problem? Fax your purchase order to 830-608-0882.

For discount on shipping multiple items call 877-608-9812.

[Link to calculate shipping for P.O.'s](#)

[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.

Prices subject to change without notice.

Turning Point Therapy & Technology,  
Inc.  
PO Box 310945  
New Braunfels, TX 78131-0945

Member of Assistive Technology Industry  
Association

### Contact Us

Toll Free: 877-608-9812  
 Phone//Fax: 830-608-9812//830-608-0882  
 Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
 Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
 Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)





Since 11/24/98



[← Back to Switches](#)


## Strong Switches

We have a wide range of single switches used for people with a wide range of capabilities, all with an 1/8 in mono plug and a 5' cable, including square pad, round pad, platform, treadle, flex, and grasp switch

Item	Description	Photo	Price
TA-58550	<p><b>Flex Switch</b></p> <ul style="list-style-type: none"> <li>• Activated by bending the black rubber tip in any direction; up, down, forward, and backward! Attached tubing is bent to facilitate positioning.</li> <li>• Tactile and auditory feedback.</li> <li>• Requires mounting.</li> </ul>		<b>\$175.00</b>
TA-58650	<p><b>Grasp Switch</b></p> <ul style="list-style-type: none"> <li>• A squeeze or pinch will activate this switch.</li> <li>• No tactile or auditory feedback.</li> <li>• The Grasp Switch is designed to be activated with a Gentle Squeeze or Pinch.</li> <li>• The user must be capable of releasing the squeeze within a short period of time.</li> <li>• If the switch is held for over 2 seconds, then air is trapped in the cylinder and the switch will latch or lock on.</li> <li>• If the user can not quickly release this pattern of movement then the Grasp Switch is not the appropriate switch for them.</li> </ul>		<b>\$137.00</b>
<p>Switch TA-58000</p> <p>Connection TA-66130</p>	<p><b>Platform Switch and Connection</b></p> <ul style="list-style-type: none"> <li>• Very durable switch which is activated by pressing on the brown raised pad.</li> <li>• Tactile and auditory feedback.</li> </ul>	<p>Switch</p>  <p>Connection</p> 	<p>Switch <b>\$99.00</b></p> <p>Connection <b>\$37.00</b></p>

TA-58100 Square Pad Switch


- Large surface area advantageous for people who have difficulty targeting a small switch.
- Activated by pressing anywhere on top surface.
- Tactile and auditory feedback.



\$99.00

TA-58800 Treadle Switch

- Great as a floor/foot switch since you can rest the activating body part on the top of the switch between activations.
- Activated by pressing down on hinged top surface.
- Tactile and auditory feedback.



\$42.00

Using PayPal a problem? Fax your purchase order to 830-608-0882.

For discount on shipping multiple items call 877-608-9812.

[Link to calculate shipping for P.O.'s](#)

[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.

Prices subject to change without notice.

Turning Point Therapy & Technology,  
Inc.  
PO Box 310945  
New Braunfels, TX 78131-0945

Member of Assistive Technology Industry  
Association

Contact Us

Toll Free: 877-608-9812  
Phone//Fax: 830-608-9812//830-608-0882  
Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)

Since 11/24/98






[•-- Back to Switches](#)

## Dual Switches

Pneumatic Switch and Rocker Switch with connector/mount and USB mouse mover

Item	Description	Photo	Price
<b>Switch</b> Short Neck TA-51000	<b>Pneumatic Switch and Mount/Connection</b> <a href="#">Link to Mounting Options</a>	Switch 	Switch <b>\$319.00</b>
Long Neck TA-51010	<ul style="list-style-type: none"> <li>- An excellent switch for people with gross motor skills.</li> <li>- A Dual Switch is two single switches in one. It has a 1/8" or 3.5mm stereo plug.</li> <li>- A sip and puff switch: to activate one switch, you "sip" and to activate the other, you "puff".</li> <li>- Internal saliva trap allows for easy cleaning.</li> <li>- <i>See Advanced Solutions below.</i></li> </ul>	Connection 	Connection <b>\$49.00</b>
Connection TA-66140	<ul style="list-style-type: none"> <li>- 2 models:                Pneumatic-short (5"/12.5cm neck)                Pneumatic-long (14"/35cm neck) with extended tube for swivel clamp mounting</li> <li>- <i>For health reasons, this product is non returnable.</i></li> </ul>		
TA-51050	<b>Rocker Switch</b> <ul style="list-style-type: none"> <li>- The switch on each side is activated by pressing down on any part of the surface.</li> <li>- <i>See Advanced Solution below.</i></li> </ul>		<b>\$131.00</b>
	<b>Advanced Solutions</b>		
TA-43410	<b>CA2 Connector</b> <p>Do you need your dual pneumatic switch to act as two single switches?            Allows you to use a dual switch with 2 devices that are each operated by a single switch</p>		<b>\$35.00</b>

TA-43540

## CA6 Connector

\$90.00

Allows you to use 1 multiple switch adapted to make 5 single switches (5 devices that are each operated by a single switch)



TA-25006

## The USB Mouse Mover

\$325.00

- A multiple switch does not plug into a computer.
- The 9 pin D female connector is not a serial connection.
- For computer access with a multiple switch, USB Mouse Mover
- Allows an individual to control all mouse functions on their computers using a combination of 6 single switches or any multiple switch.
- The single switches plug into the 1/8" jacks on the front of the Mouse Mover, which plugs directly into the computer.
- A multiple switch plugs into the 9 pin jack on the Mouse Mover.
- Pressing on each switch then delivers a specific mouse function: click, press and hold, release, and all mouse directional movements.
- Use with any software, including on-screen keyboards.



Using PayPal a problem? Fax your purchase order to 830-608-0882.

For discount on shipping multiple items call 877-608-9812.

[Link to calculate shipping for P.O.'s](#)

[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.

Prices subject to change without notice.

Turning Point Therapy & Technology,  
Inc.  
PO Box 310945  
New Braunfels, TX 78131-0945

Member of Assistive Technology Industry  
Association

## Contact Us

Toll Free: 877-608-9812  
Phone//Fax: 830-608-9812/830-608-0882  
Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)

Since 11/24/98



[← Back to Switches](#)

## Multiple Switches

Each multiple switch is 5 switches in one and comes in 1 of 2 types:

1. Those with a DB9 connector are mainly for wheelchairs; they do not function through computer serial connections without an interface, such as the [USB Mouse Mover](#)
2. Those with a USB connector function through the USB port of any computer

Item	Description	Photo	Price
DB9 Plug <a href="#">TA-50120</a>	<p><b>Joystick with Pad</b></p> <ul style="list-style-type: none"> <li>• A large, heavy duty switch for people who can maneuver a joystick but need a large one.</li> <li>• The joystick has 4 positions, and the fifth switch is activated by pressing the brown raised pad.</li> <li>• Tactile and auditory feedback.</li> <li>• See discussion at top of page regarding the 2 types of interface/plug/connector.</li> </ul> <p>Specifications:  <i>Dimensions:</i> 11 x 7 x 2.75 in (28 x 18 x 12 cm)  <i>Activating Force:</i> 6.5 oz (188 g)  <i>Travel Before Activation:</i> ball: 0.02 in (0.05 cm) / pad: 0.1 in (0.2 cm)  <i>Current:</i> 30VDC max 1A</p>	<p>DB9 Plug</p> 	<p>DB9 Plug <b>\$310.00</b></p>
USB Plug <a href="#">TA-25312</a>		<p>USB Plug</p> 	<p>USB Plug <b>\$320.00</b></p>
DB9 Plug <a href="#">TA-50200</a>	<p><b>MiniJoystick with Push</b></p> <ul style="list-style-type: none"> <li>• 4 switches are activated by moving the joystick in 4 directions; the fifth is activated by pressing down on the joystick shaft.</li> <li>• 2 tips are included: one is flexible rubber and can be cut to desired size &amp; the other is hard plastic.</li> <li>• Tactile and auditory feedback.</li> <li>• Short throws.</li> <li>• Easily mounted with the <a href="#">U clamp</a> (TA-59700) to 1/2 " tubes.</li> <li>• See discussion at top of page regarding the 2 types of interface/plug/connector.</li> </ul> <p>Specifications:  <i>Dimensions:</i> 1.25 dia x 4.5 in (3 dia x 12 cm)  <i>Activating Force:</i> joystick: 4.5 oz (128 g) / push: 7.1 oz (205 g)  <i>Travel Before Activation:</i> joystick: 0.5 in (1.3 cm) / push: 0.1 in (0.2 cm)  <i>Current:</i> 30VDC max 1A</p>	<p>DB9 Plug</p> 	<p><b>Plug Options</b>            DB9 Plug <b>\$294.00</b></p>
USB Plug <a href="#">TA-25320</a>		<p>USB Plug</p> 	<p>USB Plug <b>\$325.00</b></p> <p><b>U Clamp</b> <b>\$33.00</b></p>

DB9 Plug  
TA-50240

USB Plug  
TA-25324

### Penta Switch

- 4 directional control with center switch for mouse click.
- For people with fine motor skills and good targeting ability.
- The 5 switches are small (.3" / 1 cm diameter) and activated by pressing down on buttons.
- Tactile and auditory feedback.
- Not meant for harsh or rugged environments.
- See discussion at top of page regarding the 2 types of interface/plug/connector.

#### Specifications:

*Dimensions:* 2.0 dia x 0.5 in (5 dia x 1.3 cm)  
*Activating Force:* joystick: 4.2 oz (127 g)  
*Travel Before Activation:* 0.03 in (0.07 cm)  
*Current:* 12VDC max 50mA

DB9 Plug  
TA-50260

USB Plug  
TA-25326

### Star Switch

- Each switch is 2.25" / 6.5 cm diameter and slightly recessed.
- Comes with five colored stickers.
- Slight tactile and auditory feedback.
- See discussion at top of page regarding the 2 types of interface/plug/connector.

#### Specifications:

*Dimensions:* 9 x 9 x 0.5 in (23 x 23 x 1.3 cm)  
*Activating Force:* 5 oz (142 g)  
*Travel Before Activation:* 0.03 in (0.07 cm)  
*Current:* 12VDC max 50mA

DB9 Plug  
TA-50180

USB Plug  
TA-25318

### Wafer Switch

- 5 low profile membrane switches, slightly recessed to provide a keyguard-like effect.
- Each switch is 2.25" / 6.5 cm diameter and printed with a direction arrow
- No tactile or auditory feedback.

#### Specifications:

*Dimensions:* 15.5 x 6 x 5 in (50 x 15 x 13 cm)  
*Activating Force:* joystick: 3.5 oz (99 g)  
*Travel Before Activation:* 0.02 in (0.05 cm)  
*Current:* 24VDC max 25mA

DB9 Plug



DB9 Plug  
\$168.00

USB Plug  
\$183.00

USB Plug



DB9 Plug



DB9 Plug  
\$210.00

USB Plug  
\$225.00

USB Plug




DB9 Plug



DB9 Plug  
\$210.00

USB Plug  
\$225.00

USB Plug



Using PayPal a problem? Fax your purchase order to 830-608-0882.  
For discount on shipping multiple items call 877-608-9812.  
[Link to calculate shipping for P.O.'s](#)  
[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.  
Prices subject to change without notice.

Turning Point Therapy & Technology,  
Inc.  
PO Box 310945  
New Braunfels, TX 78131-0945

Member of Assistive Technology Industry  
Association

**Contact Us**


Toll Free: 877-608-9812  
Phone//Fax: 830-608-9812//830-608-0882  
Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)

Since 11/24/98



[← Back to Switches](#)

## Adjustable Pressure Switches

Item	Description	Photo	Price	
IT-AP	Adjustable Pressure Switches		1 Switch <b>\$59.00</b> Choose Color:	
<b>Color</b>	Some users will only be able to exert the lightest of pressure to press a switch. Others will prefer a firm pressure to minimize unwanted presses, perhaps resulting from tremor. The Adjustable Pressure Switch serves each of the users.			
Red 91				
Yellow 94				
Green 92	These robust, general purpose switches allow easy adjustment of the operating pressure. Just twist the cap and the pressure needed to activate it can be changed in small steps from light to heavy. Available in 5" diameter. The solid hexagonal bases are drilled to fit standard switch mounting plates. The back of the switch also has tappings to fit the TASH mounting system. The switch terminates in a standard 3.5mm jack plug.		All 4 Colors <b>\$199.00</b>	
Blue 93				
	Please select color when ordering. <i>Colors available:</i> red, green, blue and yellow.			

Using PayPal a problem? Fax your purchase order to 830-608-0882.

For discount on shipping multiple items call 877-608-9812.

[Link to calculate shipping for P.O.'s](#)

[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.

Prices subject to change without notice.

Turning Point Therapy & Technology,  
Inc.

PO Box 310945  
New Braunfels, TX 78131-0945

Member of Assistive Technology Industry  
Association

### Contact Us





Toll Free: 877-608-9812  
Phone//Fax: 830-608-9812/830-608-0882  
Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)

Since 11/24/98



[← Back to Mounting](#)

## Accessories, Cables, Kits

Item	Description	Photo	Price
TA-65955	<p><b>Ability Switch Tester</b></p> <ul style="list-style-type: none"> <li>• The Ability Switch Tester is a great unit for testing your switches.</li> <li>• All 3 types of switches - single, dual, and multiple - plug in directly and may be tested. • Battery included.</li> </ul> <p>Suitable for:</p> <ul style="list-style-type: none"> <li>• Evaluating a number of different clients on switches: an auditory beep and light will indicate that your client has successfully activated a switch</li> <li>• Isolating problems: if your switch will not activate a particular device, you can find out whether it's your switch or the device that is not working!</li> </ul>		<b>\$130.00</b>
CA-1 TA-43420	<p><b>CA-1 Cable Connector</b></p> <p>Allows you to plug 2 single switches into a device with a stereo jack</p>		<b>\$35.00</b>
CA-2 TA-43410	<p><b>CA-2 Cable Connector</b></p> <p>Allows you to use a dual switch with 2 devices that are each operated by a single switch</p>		<b>\$35.00</b>
CA-4 TA-43580	<p><b>CA-4 Cable Connector</b></p> <p>Allows you to use a dual switch with 1 device that is operated by a single switch</p>		<b>\$35.00</b>

CA-5  
TA-43530

### CA-5 Cable Connector

Allows you to use 5 single switches to operate any device that accepts 1 multiple switch: that has a 9 pin D sub connector



\$110.00

CA-6  
TA-43540

### CA-6 Cable Connector

Allows you to use 1 multiple switch with 5 devices that are each operated by a single switch



\$90.00

Single switch  
TA-43330

Dual switch  
TA-43350

Multiple switch  
TA-43370

### Cable Extension

- Cable extensions 5ft - 152 cm

\$25.00

TA-34500

### Call Bell

- The Call Bell is a small battery operated alarm which when activated will emit a distinctive two tone chime.
- The Call Bell can be operated by direct access, a single switch or remotely with an Ultra One package.



\$90.00

TA-65956

### Single Switch Tester

\$40.00

TA-43880

### Switch Latch & Timer

Change the action of a switch so that the device remains on even after you let go of the switch.

- Use with adapted toys, tape recorders, radios, etc.
- May use up to two switches and two separate devices.
- Modes of activation: momentary, latching, and timed.
- Battery included.



\$131.00

AA or C Size  
TA-43780

D size  
TA-43790

**Toy Cables**

- To adapt your toys for switch Access you will need:
- Any single switch ... please see single switches.
- A Toy Cable (3 feet long).
- Any battery-operated toy or appliance with AA, C, or D batteries.
- Insert the copper disk at the end of the Toy Cable between the battery and battery contact.
- Turn the toy on.
- It should not work, since you have interrupted its connection to the battery. - Plug your switch into the 1/8" jack at the other end of the Toy Cable.
- Activate the switch to operate the toy.
- Batteries sold separately.



**\$11.00**

**Kits**

TA-65990

**Deluxe Assessment Kit**  
[Link to description](#)

**\$3570.00**

TA-65992

**Mounting Assessment Kit**  
[Link to description](#)

**\$683.00**

TA-65994

**Single Switch Kit**  
[Link to description](#)

**\$1395.00**

TA-66005

**Mounting Arm Assessment Kit**  
[Link to description](#)

**\$525.00**

[Top of Page](#)

**Deluxe Assessment Kit** (TA-65990) includes the following:

Buddy Button	TA-57400
Cap	TA-58060
Cup	TA-58070
Flex	TA-58550
Grasp	TA-58650
Leaf	TA-58600
Micro Light	TA-58500
Mini Cup	TA-58080
Pillow	TA-58750
Plate-red	TA-58300
Platform	TA-58000
Ribbon	TA-58700
Round Pad	TA-58050
Soft-red	TA-58200
Square Pad	TA-58100
Tip	TA-58400
Treadle	TA-58800
Big Buddy	TA-56100
Mini Lever	TA-53100
Trigger	TA-55200
SCATIR Deluxe	TA-51300
Pneumatic-short	TA-51000

[Top of Page](#)

**Mounting Assessment Kit** (TA-65992) includes the following:

Gooseneck 13"	TA-59310
Pinch Clamp - 1"	TA-59460
Platform Clamp	TA-59490
Pinch Clamp - 7/8"	TA-59470
Switch Clamp	TA-59500
Tube - 12" w/Bushing	TA-59780
Tube - 16"	TA-59760
Tube - 30 bend, 12"	TA-59810
Tube Clamps	TA-59540
U Clamp	TA-59700
Universal Table Mount	TA-59400
Universal Swivel Clamp	TA-59450
Wheelchair Swivel Clamp	TA-59480
Mounting Arm	TA-66000
Mounting Clamp	TA-66200
Mini Cup / Trigger Connection	TA-66110
Buddy Connection	TA-66120
Big Buddy Connection	TA-66130
Switch and Mounting CD	
Carry Case	

[Top of Page](#)

[Top of Page](#)

**Single Switch Kit** (TA-65994) includes the following:

Single Switch Tester	TA-65956
Buddy Button	TA-57400
Cup	TA-58070
Grasp	TA-58650
Ribbon Switch	TA-58700
Micro Light	TA-58500
Mini Cup	TA-58080
Pillow	TA-58750
Plate-red	TA-58300
Soft-green	TA-58210
Square Pad	TA-58100
Switch Latch	TA-43880
Tip	TA-58400
Treadle	TA-58800
Cap	TA-58060
Toy Cable AA/C Battery	TA-43780
Toy Cable D Battery	TA-43790
Big Buddy	TA-56100
Trigger	TA-55400
Mini Lever	TA-53110
SCATIR Deluxe	TA-51300
Switch and Mounting CD	

Rocker	TA-51050
Joystick with Pad	TA-50120
Mini Joystick with Push	TA-50200
Penta	TA-50240
Star	TA-50260
Wafer	TA-50180
Ability Switch Tester	TA-65955
CA-1	TA-43420
CA-2	TA-43410
CA-4	TA-43580
CA-5	TA-43530
CA-6	TA-43540
Switch Latch	TA-43880
Toy Cable AA/C battery	TA-43780
Toy Cable D battery	TA-43790
Mounting Arm	TA-66000
Mounting Clamp	TA-66200
Mini / Trigger Connection	TA-66110
Buddy Connection	TA-66120
Ribbon Connection	TA-66150
Big Buddy Connection	TA-66130
Microlight Connection	TA-66100
Switch and Mounting CD	
Carry Case	

**Mounting Arm Assessment Kit** (TA-66005) includes the following:

Mounting Arm	TA-66000
Microlight Connection	TA-66100
Mini Connection	TA-66110
Buddy Connection	TA-66120
Big Buddy Connection	TA-66130
Pneumatic Connection	TA-66140
Ribbon Connection	TA-66150
Relax Mini Relax Connection	TA-66160
Mounting Clamp	TA-66200
Universal Clamp Connection	TA-66210
1/2 in Bottom Tube Connection	TA-66230
1/2 in Top Tube Connection	TA-66240
Switch and Mounting CD	
Carry Case	

Carry Case

Using PayPal a problem? Fax your purchase order to 830-608-0882.

For discount on shipping multiple items call 877-608-9812.

[Link to calculate shipping for P.O.'s](#)

[Shipping outside the continental United States](#) requires [email contact](#) prior to purchase.

Prices subject to change without notice.

Turning Point Therapy & Technology, Inc.  
PO Box 310945  
New Braunfels, TX 78131-0945

Member of Assistive Technology Industry Association

#### Contact Us

Toll Free: [877-608-9812](tel:877-608-9812)  
Phone//Fax: [830-608-9812](tel:830-608-9812)//[830-608-0882](tel:830-608-0882)  
Support: [support@turningpointtechnology.com](mailto:support@turningpointtechnology.com)  
Sales: [sales@turningpointtechnology.com](mailto:sales@turningpointtechnology.com)  
Website: [webmaster@turningpointtechnology.com](mailto:webmaster@turningpointtechnology.com)

Since 11/24/98

