




Communication

Speech output & scanning devices & picture symbol board

[← Back to Switches, Computer Access, Controls](#)

Item	Description	Photo	Price
IT-FT04	<p>Traxsys 4Talk4</p> <p>A durable classroom speech output device from Traxsys.</p> <p>4 large colored buttons, each with its own message. 4 distinct levels for pre-recorded messages. Plastic pocket inserts allow you to use your own pictures to represent the messages.</p> <p>Clear, high-quality speech with adjustable volume control (loud enough for school assembly!). In addition, the 4talk4 has ports to allow you to connect external switches for those who need to use their own switches, and ports so you can attach and control battery operated toys and devices. Press to say "Make the rabbit jump" and then see it happen! Switches and toys not included but the batteries are.</p>		\$349.00
IT-SC71	<p>Traxsys Scan4</p> <p>Scan4 is a scanning version of the popular 4talk4 which has been carefully designed to be used as an introduction to scanning. Scan4 has very bright scanning lights, auditory scanning and can use single or two switch scanning. As with 4talk4, it has sockets to connect battery operated devices, and comes with two toy leads.</p> <p>Encourage your switch users to understand scanning and selection of their choice. Bright lights, clear speech and a real outcome (with toys connected).</p> <p>Traditional single switch scanning, starting at a realistically slow speed. Two switch scanning for the user who can manage two switches and have more control.</p> <p>It also includes a single switch "dwell" scan mode, where pressing the switch controls the scan, and dwelling on the one you want causes it to activate.</p> <p>In direct mode, press the area for the message. The light above the area shines to reinforce which one has been pressed.</p>		\$599.00
ComBoard TA-26150	<p>ComBoard</p> <ul style="list-style-type: none">• The Comboard is a picture symbol board which houses a motor-driven pointer.• The pointer is rotated until it reaches the desired picture symbol on the board by pressing on 1 of 2 single switches which plug directly into 1/8" (3.5 mm) jacks on the side of the unit.• One switch rotates the pointer in one direction, the other in the opposite direction.• The speed at which the pointer rotates can be adjusted.• Suitable for beginning communicators who require special access using single switches• <i>Single Switches and Extra Picture Symbols Sold Separately</i>		ComBoard \$250.00
Extra Symbols TA-26070			Extra Symbols \$10.00

[Details](#)

Red
TA-27650

Yellow
TA-27660

Talking Buddy

- 20 seconds of recorded time
- 9V battery included
- Option to use alternate external single switch
- Volume on-off control Cable connection jack allows you to access any switch adapted toy or appliance at the same time
- Available in **yellow** (shown to right) or **red**

Specifications:

Size: 4.5 x 1.5" / 1 lb

Weight: 11 x 4 cm / 300 g



\$110.00

Call for discount on shipping multiple items.

Shipping outside the United States requires email contact prior to purchase.

Prices subject to change without notice.

Turning Point Therapy & Technology, Inc.

PO Box 310945

New Braunfels, TX 78131-0945

Contact Us

Toll Free: 877-608-9812

Phone/Fax: 830-608-9812//830-608-0882

Support: support@turningpointtechnology.com

Sales: sales@turningpointtechnology.com

Website: webmaster@turningpointtechnology.com

Member of Assistive Technology Industry Association



Since 11/24/98

